

# Tap Water Iontophoresis Device SweatStop® DE20



## Instruction for Use Head Electrode

Before putting on the head electrode, moisten it strongly, to the extent that it barely does not drip. Then, connect the electrode cable to the two electrodes (black pads). Put on the head electrode like a hood and fasten all Velcro straps tightly.



Electrode

**It is important to ensure that the electrodes on the neck and forehead are very flat on the skin and that there are not only punctual contact points. If necessary, tighten the Velcro straps or realign the hood.**

Before turning on the device, connect the cables to the device. For the further steps, please refer to the instruction manual enclosed with the device.

**The current for the treatment of the head should not be set higher than 5 mA. The ideal intensity of current for the facial treatment lies in between 1 to 3 mA (approximately level 0–1.5 for GS, or level 2–4 for PS).**

If it should medically be necessary to set the current higher than 2 mA, this must be done under the supervision of a doctor and possible skin irritations are to be carefully monitored.

As a rule of thumb for the optimal amperage, turn up the current until you can feel the current, then turn it back a little so that you can no longer feel the current flow (GS, direct current) or slightly feel it (PS, pulsating current).

### Indication

If the current cannot be activated, the resistance between the electrodes is too high and the protective circuit against over-therapy is activated. In this case you need to remove the facial electrode again. Fill a tub with lukewarm water and add a teaspoon of table salt to it.

Stir the water until the salt is completely dissolved. Moisten the head electrode in this tub and apply it anew. Start the treatment again as described above.